


ADP Glo™ assay

Liesbeth Steegmans, ing. Cistim, Leuven, Belgium



Introduction

ADP Glo™ assay

We use an ADP-Glo™ assay (Promega) (see **Figure 1**) to test different compounds in dose response and in the meantime look at the substrate dependency of this compound. We work in a relative small assay volume so to make the pipetting more doable it is easy to spot a dose response of the substrate in de assay plate, and add the compound dose response manually or with a robotic system (Fluent).

Materials and methods

We make a dose response of the substrate in a buffer solution. We place the cups in the source plate of the I.DOT and fill them with the dose response of the substrate. We spot this dose response with the I.DOT on the destination plate (see **Figure 2**). Then the compound dose response is added manually or with a robotic system like the Fluent (TECAN). The reaction is started by dropping ATP using a multi-drop, this allows us to work fast, robust and gives the opportunity to drop small volumes.




Figure 1. ADP Glo™ assay principle

APPLICATION NOTE

After 30 min at 37°C the reaction is stopped with ADP Glo™ Stop using a multidrop to deplete not consumed ATP. After 40 min at RT ADP Glo™ Detection is added with a multidrop to the plate to

convert ADP to ATP. After 1h at RT we measure Luminescence with an Envision. All the data is plotted with Prism (GraphPad)




Figure 2. Layout from the software

Results

Graphs of the Dose dependency of the Substrate and the shift (Vmax/Km) with different compound concentrations is calculated (see Figure 3).

Conclusions

With the I.DOT it is possible to spot easy and precise the dose response of the substrate. So that we have more volume left to add other components to the reaction. This set up gives us the possibility to get the most out of this experiment. So we can decide how to proceed with this project.




Figure 3. Graphs from graphpad

References

1. https://be.promega.com/-/media/files/resources/protocols/technical-manuals/0/adp-glo-kinase-assay-protocol.pdf?rev=42e85ccab87244b395f4ecd1ad8e1616&sc_lang=en



©2021 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.