


I.DOT & DropDetection: A unique feature for quality control and dispensing protocol optimization.

Alexander Jurkowitsch & Trang Nguyen, DISPENDIX GmbH, Stuttgart, Germany



Introduction

Non-contact, low-volume dispensing depends on generating individual and discrete droplets. They typically range from pico- to microliters, and volumes this small can't be evaluated visually. Immediate Drop-On-Demand Technology (I.DOT) is a non-contact, low-volume dispensing platform. The general principle is based on a hole (less than 100 μm) in the bottom of a microtiter plate well, where surface tension and capillary forces retain the sample liquid in the cavity. Applying a pulse of air pressure on top of the well forms a precise droplet that can be released into or onto almost any target.

Each I.DOT is equipped with DropDetection to interface with the low-volume world.

DropDetection is a patented feature unique to the I.DOT. It enables detection and counting of every droplet released during a dispensing run. It acts as a simple yet powerful tool for quality control and protocol optimization.

DropDetection is mounted under the I.DOT One source tray (**Figure 1**). It's a circuit board that leverages 96 miniaturized light barriers to detect every droplet that each source plate position generates. DropDetection identifies changes in light intensity to detect droplets passing the light barrier. After dispensing is complete, DropDetection provides a text file and color-coded result. Green indicates successful dispensing and red indicates an error. If the software recognizes the error, it displays the potential reason as well.




Figure 1. Individual droplets dispensed through source tray of I.DOT One and DropDetection (1) circuit board underneath source tray of I.DOT One, LED (2), photodiode (3) and dispensed liquid (4).

APPLICATION NOTE

DropDetection output notes:

- Passed – everything went accordingly.
- Strange – unexplainable error.
- Unfilled/Clogged – the well wasn't filled to begin with, or got clogged.
- Became empty – the liquid in the well became empty; there was not enough for the completion of the assay.

We designed an experiment to validate DropDetection's performance.

Experiment design

We setup a dispensing protocol to demonstrate four scenarios (**Figure 2**). Four liquids were dispensed into a 384-well target plate.

"Water 1" (**yellow**) represented normal water.

Fifty μL were added to the source well. No error was expected.

"Water 2" (**green**) also represented water. Here, an insufficient amount (less than $14.4 \mu\text{L}$ as calculated by the software) to run the assay was added to the source well. The error "became empty" was expected.

"Water 7" (**blue**) represents water mixed with a foaming agent to mimic unreliable dispensing behavior. The error "strange" was expected.

"Water 5" (**purple**) is not filled. The error "unfilled/clogged" was expected.

After setting up and filling the source wells with liquids, the I.DOT was ready for dispensing. To dispense a total volume of 150 nL per target well, 3 droplets of 50 nL droplet volume are generated.




Figure 2. Example setup of a dispensing protocol. (1) source plate layout, (2) liquids used in the setup (with indication of minimum volume to complete the dispensing task), (3) target plate layout.




Figure 3. Color-coded DropDetection results.

Results

The results demonstrate the I.DOT's dispensing accuracy in each scenario (**Figure 3**).

"Water 1" was dispensed successfully. The software indicates this with the performance note "passed."

Three droplets of 50 nL were expected and dispensed in well B20.

"Water 2", as expected, began dispensing successfully. Due to insufficient volume in the dispensing well, the well ran empty. The software indicates this with the performance note "became empty." Three droplets of 50 nL were expected in well H16, but only 1 droplet was dispensed.

For "Water 7", the software shows the performance note "strange," and the color scheme notifies the user of a dispensing error. Three droplets of 50 nL were expected in well J10, but only 2 were dispensed, possibly due to air bubbles produced by the foaming agent.

For "Water 5", the software indicates that the source well was "unfilled" or "clogged" from the beginning. Three 50 nL droplets were expected and none were released. Each dispensing run generates a DropDetection report file (**Figure 4**). The report log displays the dispensed liquid, positions of source and target wells and the detected droplets of each individual liquid-transfer step. These report files can be used to evaluate dispensing performance, and act as immediate quality control.

APPLICATION NOTE

```
DateTime: 2019-04-10 14:46:45
SoftwareVersion: 1.3.20180912

Source Plate ID: Source Plate 1 (PurePlate)
Target Plate ID: Target Plate 1 (MMP 384)

Liquid      Source Target Drop Miss TargetVolume
Water 1     A2    A1    3    0    150 nL
Water 2     C2    E1    3    0    150 nL
Water 7     E2    I1    3    0    150 nL
Water 5     G2    M1    3    0    150 nL
Water 1     A2    B1    3    0    150 nL
Water 2     C2    F1    3    0    150 nL
Water 7     E2    J1    3    0    150 nL
Water 5     G2    N1    3    0    150 nL
Water 1     A2    C1    3    0    150 nL
Water 2     C2    G1    3    0    150 nL
Water 7     E2    K1    3    0    150 nL
Water 5     G2    O1    3    0    150 nL
Water 1     A2    D1    3    0    150 nL
Water 2     C2    H1    3    0    150 nL
Water 7     E2    L1    3    0    150 nL
Water 5     G2    P1    3    0    150 nL
Water 1     A2    A2    3    0    150 nL
Water 2     C2    E2    3    0    150 nL
Water 7     E2    I2    3    0    150 nL
```

Figure 4. DropDetection report file

Conclusion

DropDetection reliably and robustly detected the expected outcome.

This example was performed with water. However, the DropDetection works with every dispensable liquid — including DMSO, glycerol-containing solutions

and more. Contact Dispindex for more information on dispensable liquids.

DropDetection has great potential to control dispensing quality. It's a noninvasive recording technology that accurately identifies dispensing faults to optimize dispensing protocols.



©2021 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.